



Rabbit Anti-Arabidopsis thaliana CUL3A Polyclonal Antibody (CABT-Z977R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH-conjugated synthetic peptide of CUL3A derived from Arabidopsis thaliana AT1G26830.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Unpurified
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1:1000-1:2000
Molecular Weight	85 kDa
Reconstitution	Reconstitution with 150 µl of sterile water.
Format	Lyophilized
Size	150 µg
Preservative	None
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 6 months, -20 to -70 °C under sterile conditions after reconstitution. 1 month, 2 to 8 °C under sterile conditions after reconstitution.
Ship	Wet ice

BACKGROUND

Introduction	CUL3A is ubiquitously expressed in plants and is able to interact with the ring-finger protein RBX1. Arabidopsis CUL3 forms E3 protein complexes with certain BTB domain proteins.
Keywords	CUL3A;ATCUL3;ATCUL3A;CUL3;CULLIN 3;CULLIN 3A.

GENE INFORMATION

Gene Name	CUL3A
Entrez Gene ID	839226
UniProt ID	Q9ZVH4