



Human Anti-Human beta Amyloid 42 (aa 1-5) monoclonal antibody, Clone A6 (CABT-ZB006)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Bind to human A β
Target	beta Amyloid
Isotype	IgG
Source/Host	Human
Species Reactivity	Human
Clone	A6
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA
Epitope	aa 1-5
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Beta-amyloid is a 40–43 amino acid peptide cleaved from amyloid precursor protein, and a putative gamma-secretase. It is a major component of the extracellular plaques found in Alzheimer's disease brain tissue.
Keywords	A4_HUMAN; AAA; ABETA; ABPP; AD1; AICD-51

GENE INFORMATION

Synonyms	Aβ; Abeta
Entrez Gene ID	351
UniProt ID	P05067