



# Mouse Anti-Clobetasol propionate monoclonal antibody, clone CBP (CABT-ZB012)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Clobetasol propionate
Immunogen	Clobetasol propionate with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	CBP
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

# BACKGROUND

Introduction	Clobetasol propionate is the 17-O-propionate ester of clobetasol. A potent corticosteroid, it is used to treat various skin disorders, including exzema and psoriasis. It has a role as an anti-inflammatory drug. It is an 11beta-hydroxy steroid, a 20-oxo steroid, a glucocorticoid, a fluorinated steroid, a 3-oxo-Delta(1),Delta(4)-steroid and a chlorinated steroid. It derives from a clobetasol and a propionic acid.
--------------	--

Keywords	CLOBETASOL PROPIONATE; 25122-46-7; Clobetasol 17-propionate; Clobex; Dermovate
----------	--

# GENE INFORMATION

Synonyms	CLOBETASOL PROPIONATE; 25122-46-7; Clobetasol 17-propionate; Clobex; Dermovate
----------	--