



Mouse Anti-11 α -hydroxyprogesterone monoclonal antibody, clone 11AOHP (CABT-ZB013)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	No cross-reactivity with progesterone and 17 α -hydroxyprogesterone
Target	Hydroxyprogesterone
Immunogen	11 α -hydroxyprogesterone with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	11AOHP
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None

Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Hydroxyprogesterone is a form of progestin, a manmade form of a female hormone called progesterone.
Keywords	11-Hydroxyprogesterone 11-glucuronide; (2S,3S,4S,5R,6R)-6-[[[(10R,11R,13S,17S)-17-acetyl-10,13-dimethyl-3-oxo-1,2,6,7,8,9,11,12,14,15,16,17-dodecahydrocyclopenta[a]phenanthren-11-yl]oxy]-3,4,5-trihydroxyoxane-2-carboxylic acid; 11-Hpg; 129238-72-8; 11alpha-Glucuronylprogesterone

GENE INFORMATION

Synonyms	11-Hydroxyprogesterone 11-glucuronide; (2S,3S,4S,5R,6R)-6-[[[(10R,11R,13S,17S)-17-acetyl-10,13-dimethyl-3-oxo-1,2,6,7,8,9,11,12,14,15,16,17-dodecahydrocyclopenta[a]phenanthren-11-yl]oxy]-3,4,5-trihydroxyoxane-2-carboxylic acid; 11-Hpg; 129238-72-8; 11alpha-Glucuronylprogesterone
----------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------