



# Mouse Anti-Acid scarlet G monoclonal antibody, clone ASG (CABT-ZB020)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	Acid red 97
<b>Immunogen</b>	Acid scarlet G with carrier protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	ASG
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, LFIA
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	Long time storage is recommended at -20°C.
<b>Ship</b>	Wet ice

# BACKGROUND

**Introduction** C.i. acid red 97 is a red to dark red powder. Copper and iron dull color; iron to a greater extent.

---

**Keywords** Acid scarlet G; ASG, Acid Red 97; 10169-02-5; C.I. Acid Red 97; UNII-8UJ772G61O

---

# GENE INFORMATION

**Synonyms** Acid scarlet G; ASG, Acid Red 97; 10169-02-5; C.I. Acid Red 97; UNII-8UJ772G61O

---