



Mouse Anti-Profenofos monoclonal antibody, clone PRO (CABT-ZB022)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Profenofos
Immunogen	Profenofos with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	PRO
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Profenofos is an organic thiophosphate, an organophosphate insecticide, an organochlorine insecticide and a member of monochlorobenzenes. It has a role as an EC 3.1.1.7 (acetylcholinesterase) inhibitor, an acaricide and an agrochemical. It derives from a 4-bromo-2-chlorophenol.

Keywords

PROFENOFOS; 41198-08-7; Curacron; Selecron; Polycron

GENE INFORMATION

Synonyms

PROFENOFOS; 41198-08-7; Curacron; Selecron; Polycron