



Mouse Anti-Etoxazole monoclonal antibody, clone ETX (CABT-ZB023)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Ettoxazole
Immunogen	Ettoxazole with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	ETX
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Ettoxazole is an organofluorine acaricide.
Keywords	Ettoxazole,153233-91-1; 4-(4-(tert-Butyl)-2-ethoxyphenyl)-2-(2,6-difluorophenyl)-4,5-dihydrooxazole; 4-(4-tert-butyl-2-ethoxyphenyl)-2-(2,6-difluorophenyl)-4,5-dihydro-1,3-oxazole; 2-(2,6-Difluorophenyl)-4-(4-(1,1-dimethylethyl)-2-ethoxyphenyl)-4,5-dihydrooxazole

GENE INFORMATION

Synonyms	Ettoxazole,153233-91-1; 4-(4-(tert-Butyl)-2-ethoxyphenyl)-2-(2,6-difluorophenyl)-4,5-dihydrooxazole; 4-(4-tert-butyl-2-ethoxyphenyl)-2-(2,6-difluorophenyl)-4,5-dihydro-1,3-oxazole; 2-(2,6-Difluorophenyl)-4-(4-(1,1-dimethylethyl)-2-ethoxyphenyl)-4,5-dihydrooxazole
----------	---