



Mouse Anti-Methomyl monoclonal antibody, clone MTM (CABT-ZB024)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Methomyl
Immunogen	Methomyl with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	MTM
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Methomyl is a white crystalline solid with slight sulfurous smell. Used as a nematocide, and an insecticide on vegetables, tobacco, cotton, alfalfa, soy beans, and corn.
Keywords	methomyl; Lannate; Mesomile; Nudrin; Methomex

GENE INFORMATION

Synonyms	methomyl; Lannate; Mesomile; Nudrin; Methomex
----------	---