



Mouse Anti-Prometryn monoclonal antibody, clone PMT (CABT-ZB026)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Prometryn
Immunogen	Prometryn with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	PMT
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Prometryn is a diamino-1,3,5-triazine that is N,N'-di(propan-2-yl)-1,3,5-triazine-2,4-diamine substituted by a methylsulfanediyl group at position 6. It has a role as a herbicide, a xenobiotic and an environmental contaminant. It is a diamino-1,3,5-triazine and a methylthio-1,3,5-triazine.
Keywords	Prometryn; 7287-19-6; Prometryne; Gesagard; Prometrex

GENE INFORMATION

Synonyms	Prometryn; 7287-19-6; Prometryne; Gesagard; Prometrex
----------	---