



# Mouse Anti-1-Naphthaleneacetic acid monoclonal antibody, clone NAA (CABT- ZB027)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	NAA
Immunogen	1-Naphthaleneacetic acid with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	NAA
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.

# BACKGROUND

Introduction	1-naphthaleneacetic acid is a naphthylacetic acid substituted by a carboxymethyl group at position 1. It has a role as a synthetic auxin. It is a conjugate acid of a 1-naphthaleneacetate.
Keywords	1-Naphthylacetic acid; 86-87-3; 1-NAPHTHALENEACETIC ACID; 2-(Naphthalen-1-Yl)Acetic Acid; Naphthalene-1-acetic acid

# GENE INFORMATION

Synonyms	1-Naphthylacetic acid; 86-87-3; 1-NAPHTHALENEACETIC ACID; 2-(Naphthalen-1-Yl)Acetic Acid; Naphthalene-1-acetic acid
----------	---------------------------------------------------------------------------------------------------------------------