



# Mouse Anti-Maleic hydrazide monoclonal antibody, clone MH (CABT-ZB028)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	MH
Immunogen	Maleic hydrazide with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	MH
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

# BACKGROUND

Introduction	Maleic hydrazide is a pyridazinone. Maleic hydrazide is an odorless white solid. Sinks in water.
Keywords	MALEIC HYDRAZIDE; Pyridazine-3,6-diol; 123-33-1; 3,6-Dihydroxypyridazine; 3,6-Pyridazinediol

# GENE INFORMATION

Synonyms	MALEIC HYDRAZIDE; Pyridazine-3,6-diol; 123-33-1; 3,6-Dihydroxypyridazine; 3,6-Pyridazinediol
----------	--