



Mouse Anti-Abscisic acid monoclonal antibody, clone ABA (CABT-ZB029)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|--|
| Target | ABA |
| Immunogen | Abscisic acid with carrier protein |
| Isotype | IgG |
| Source/Host | Mouse |
| Species Reactivity | N/A |
| Clone | ABA |
| Purification | Purified from mouse ascites. |
| Conjugate | Unconjugated |
| Applications | ELISA, LFIA |
| Format | Purified, Liquid |
| Concentration | Lot specific |
| Size | 1 mg |
| Buffer | PBS |
| Preservative | None |
| Storage | Long time storage is recommended at -20°C. |
| Ship | Wet ice |

BACKGROUND

| | |
|--------------|--|
| Introduction | Abscisic acid (ABA) is a plant hormone that regulates numerous aspects of plant growth, development, and stress responses. |
| Keywords | Dormin; 14375-45-2; (+/-)-Absciscic acid; Abscisinsaeure; (RS)-Absciscic acid |

GENE INFORMATION

| | |
|----------|---|
| Synonyms | Dormin; 14375-45-2; (+/-)-Absciscic acid; Abscisinsaeure; (RS)-Absciscic acid |
|----------|---|