



Mouse Anti-Zoxamide monoclonal antibody, clone ZXM (CABT-ZB030)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Zoxamide
Immunogen	Zoxamide with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	ZXM
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Zoxamide is a member of the class of benzamides obtained by formal condensation of the carboxy group of 3,5-dichloro-4-methylbenzamide with the amino group of 3-amino-1-chloro-3-methylpentan-2-one. It is a member of benzamides, an alpha-chloroketone and a dichlorobenzene.

Keywords

156052-68-5; Zoxium; RH 7281; 3,5-dichloro-N-(1-chloro-3-methyl-2-oxopentan-3-yl)-4-methylbenzamide

GENE INFORMATION

Synonyms

156052-68-5; Zoxium; RH 7281; 3,5-dichloro-N-(1-chloro-3-methyl-2-oxopentan-3-yl)-4-methylbenzamide