



# Mouse Anti-Flubendiamide monoclonal antibody, clone FBD (CABT-ZB033)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                            |
|--------------------|--------------------------------------------|
| Target             | Flubendiamide                              |
| Immunogen          | Flubendiamide with carrier protein         |
| Isotype            | IgG                                        |
| Source/Host        | Mouse                                      |
| Species Reactivity | N/A                                        |
| Clone              | FBD                                        |
| Purification       | Affinity purification                      |
| Conjugate          | Unconjugated                               |
| Applications       | ELISA, LFIA                                |
| Format             | Purified, Liquid                           |
| Concentration      | Lot specific                               |
| Size               | 1 mg                                       |
| Buffer             | PBS                                        |
| Preservative       | None                                       |
| Storage            | Long time storage is recommended at -20°C. |
| Ship               | Wet ice                                    |

# BACKGROUND

|              |                                                                                                                                                                                                                                                                                                   |
|--------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Introduction | Flubendiamide is an organofluorine insecticide. It has a role as a ryanodine receptor modulator. It derives from a phthalamide.                                                                                                                                                                   |
| Keywords     | Flubendiamide; 272451-65-7; 3-Iodo-N2-(2-methyl-1-(methylsulfonyl)propan-2-yl)-N1-(2-methyl-4-(perfluoropropan-2-yl)phenyl)phthalamide; UNII-GEV84ZI4K6; 1-N-[4-(1,1,1,2,3,3,3-heptafluoropropan-2-yl)-2-methylphenyl]-3-iodo-2-N-(2-methyl-1-methylsulfonylpropan-2-yl)benzene-1,2-dicarboxamide |

# GENE INFORMATION

|          |                                                                                                                                                                                                                                                                                                   |
|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Synonyms | Flubendiamide; 272451-65-7; 3-Iodo-N2-(2-methyl-1-(methylsulfonyl)propan-2-yl)-N1-(2-methyl-4-(perfluoropropan-2-yl)phenyl)phthalamide; UNII-GEV84ZI4K6; 1-N-[4-(1,1,1,2,3,3,3-heptafluoropropan-2-yl)-2-methylphenyl]-3-iodo-2-N-(2-methyl-1-methylsulfonylpropan-2-yl)benzene-1,2-dicarboxamide |
|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|