



Mouse Anti-Diamide monoclonal antibody, clone DM II (CABT-ZB035)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross-Reactivity: Cyclaniliprole 0.5 ppb, chlorantraniliprole 0.5ppb, Tetrachlorantraniliprole 0.5ppb
Target	Diamide
Immunogen	Diamide with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	DM II
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None

Storage	Long time storage is recommended at -20°C.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	A sulfhydryl reagent which oxidizes sulfhydryl groups to the disulfide form. It is a radiation-sensitizing agent of anoxic bacterial and mammalian cells.
---------------------	---

Keywords	DIAMIDE; N,N,N',N'-Tetramethylazodicarboxamide; 1,1'-Azobis(N,N-dimethylformamide); 10465-78-8; TMAD
-----------------	--

GENE INFORMATION

Synonyms	DIAMIDE; N,N,N',N'-Tetramethylazodicarboxamide; 1,1'-Azobis(N,N-dimethylformamide); 10465-78-8; TMAD
-----------------	--