



# Rabbit Anti-Human DNASE1 (aa19-259) polyclonal antibody [Biotin] (CABT-ZB040)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	It has been selected for its ability to recognize DNASE1 in immunohistochemical staining and western blotting.
<b>Target</b>	DNASE1
<b>Immunogen</b>	Recombinant DNASE1 (Gly19-Ala259) expressed in E. coli.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Biotin
<b>Applications</b>	Recommended Dilutions: Western blot (1:100 - 1:400) Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	100 µL
<b>Buffer</b>	PBS with 50% Glycerol

<b>Preservative</b>	None
<b>Storage</b>	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles. Store undiluted.
<b>Ship</b>	Wet ice

## GENE INFORMATION

<b>Synonyms</b>	DNASE1; deoxyribonuclease I; DNL1; DRNI; deoxyribonuclease-1; Dornase alfa; DNase I, lysosomal; human urine deoxyribonuclease I;
<b>Entrez Gene ID</b>	<a href="#">1773</a>
<b>UniProt ID</b>	<a href="#">P24855</a>