



Rabbit Anti-Human DNASE1 polyclonal antibody (CABT-ZB041)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Predicted cross-reactivity: Rabbit, Dog(>80% identity)
Target	DNASE1
Immunogen	Recombinant protein corresponding to human DNase I
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Monkey
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	<p>WB, ICC/IF, IHC-P</p> <p>Recommended Dilutions:</p> <p>WB 1:500-1:3000</p> <p>ICC/IF 1:100-1:1000</p> <p>IHC-P 1:100-1:1000</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µL

Buffer	PBS with 20% Glycerol
Preservative	None
Storage	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles. Store undiluted.
Ship	Wet ice

GENE INFORMATION

Synonyms	DNASE1; deoxyribonuclease I; DNL1; DRNI; deoxyribonuclease-1; Dornase alfa; DNase I, lysosomal; human urine deoxyribonuclease I;
Entrez Gene ID	1773
UniProt ID	P24855