



Rabbit Anti-Human DNASE1 polyclonal antibody (CABT-ZB045)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	DNASE1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µL
Buffer	PBS with 50% Glycerol
Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

GENE INFORMATION

Synonyms	DNASE1; deoxyribonuclease I; DNL1; DRNI; deoxyribonuclease-1; Dornase alfa; DNase I, lysosomal; human urine deoxyribonuclease I;
-----------------	--

Entrez Gene ID [1773](#)

UniProt ID [P24855](#)
