



Rabbit Anti-Human DNASE1 polyclonal antibody (CABT-ZB046)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Specificity | Recognizes endogenous levels of human DNASE1. Species Crossreactivity: mouse, rat |
| Target | DNASE1 |
| Immunogen | Recombinant full length Human DNASE1. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | IHC, WB Recommended Dilutions: Western Blot: 1:500-1:2000 Immunofluorescence/Immunohistochemistry: 1:50-1:200 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. |
| Format | Purified, Liquid |
| Concentration | Lot specific |
| Size | 100 µL |
| Buffer | PBS with 50% Glycerol |

| | |
|---------------------|---|
| Preservative | 0.09% sodium azide |
| Storage | Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles. Store undiluted. |
| Ship | Wet ice |

GENE INFORMATION

| | |
|-----------------------|--|
| Synonyms | DNASE1; deoxyribonuclease I; DNL1; DRNI; deoxyribonuclease-1; Dornase alfa; DNase I, lysosomal; human urine deoxyribonuclease I; |
| Entrez Gene ID | 1773 |
| UniProt ID | P24855 |