



Rabbit Anti-Human DNASE1 polyclonal antibody (CABT-ZB047)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Immunogen Sequence Homology: Human: 100%; Rat: 79%; Bovine: 79%
Target	DNASE1
Immunogen	A synthetic peptide directed towards the N terminal of human DNASE1. Synthetic peptide located within the following region: GKLLDNLNQDAPDTYHYVVSEPLGRNSYKERYLFVYRPDQVSAVDSYYD
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µL
Buffer	PBS with 2% sucrose
Preservative	0.09% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

GENE INFORMATION

Synonyms	DNASE1; deoxyribonuclease I; DNL1; DRNI; deoxyribonuclease-1; Dornase alfa; DNase I, lysosomal; human urine deoxyribonuclease I;
Entrez Gene ID	1773
UniProt ID	P24855