



Mouse Anti-Metsulfuron monoclonal antibody, clone MSF (CABT-ZB056)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Metsulfuron
Immunogen	Metsulfuron with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	MSF
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Metsulfuron is known fully as metsulfuron methyl but is more commonly known by the abbreviations MSM. It is from the sulfonyleurea family of herbicides and is an organic active ingredient that can deliver pre-emergent and post-emergent control of broadleaf weeds and a number of annual grasses.
Keywords	Metsulfuron; 79510-48-8; UNII-2G9B203G2P; 2-[(4-methoxy-6-methyl-1,3,5-triazin-2-yl)carbamoylsulfamoyl]benzoic acid

GENE INFORMATION

Synonyms	Metsulfuron; 79510-48-8; UNII-2G9B203G2P; 2-[(4-methoxy-6-methyl-1,3,5-triazin-2-yl)carbamoylsulfamoyl]benzoic acid
----------	---