



Rabbit Anti-Human Pyrophosphatase polyclonal antibody (CABT-ZB061)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Predicted cross-reactivity: Mouse, Rat, Bovine(>80% identity)
Target	PPA1
Immunogen	Recombinant protein encompassing a sequence within the center region of human PPA1. The exact sequence is proprietary.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	Recommended Dilutions: Western blot (1:5000-1:20000) Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µL
Buffer	0.1 M Tris-Glycine (pH7) with 20% Glycerol

Preservative	None
Storage	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles. Store undiluted.
Ship	Wet ice

GENE INFORMATION

Synonyms	PPA1; pyrophosphatase (inorganic) 1; PP; PP1; IOPPP; HEL-S-66p; SID6-8061; inorganic pyrophosphatase; PPase; pyrophosphatase 2
Entrez Gene ID	5464
UniProt ID	Q15181