



Rabbit Anti-Human Pyrophosphatase (aa 1-289) polyclonal antibody [Biotin] (CABT-ZB062)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	PPA1
Immunogen	Human PPA1 (1-289 aa) expressed in E. coli
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purification
Conjugate	Biotin
Applications	ELISA, RIA
Format	Purified, Liquid
Concentration	Lot specific
Size	150 µg
Buffer	PBS with 50% Glycerol
Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

GENE INFORMATION

Synonyms	PPA1; pyrophosphatase (inorganic) 1; PP; PP1; IOPPP; HEL-S-66p; SID6-8061; inorganic pyrophosphatase; PPase; pyrophosphatase 3
Entrez Gene ID	5464
UniProt ID	Q15181