



Mouse Anti-Dieldrin monoclonal antibody, clone DLD (CABT-ZB065)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross-Reactivity: Aldrin 50 ppb
Target	Dieldrin
Immunogen	Dieldrin with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	DLD
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.

BACKGROUND

Introduction

Dieldrin is a light-tan flaked solid. It is insoluble in water. It is toxic by inhalation, skin absorption and ingestion. It can penetrate intact skin. It is used as an insecticide.

Keywords

Dieldrin; 60-57-1; HEOD; Compound 497; Octalox

GENE INFORMATION

Synonyms

Dieldrin; 60-57-1; HEOD; Compound 497; Octalox
