



Mouse Anti-Diclofop-methyl monoclonal antibody, clone DM (CABT-ZB069)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	DM
Immunogen	Diclofop-methyl with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	DM
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Diclofop methyl appears as colorless crystals. Decomposed by either strong acid or base. Used as a selective herbicide.
Keywords	Diclofop methyl; 51338-27-3; Illoxan; Dichlorfop-methyl

GENE INFORMATION

Synonyms	Diclofop methyl; 51338-27-3; Illoxan; Dichlorfop-methyl
----------	---