



Mouse Anti-Lactofen monoclonal antibody, clone LTF (CABT-ZB070)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Lactofen
Immunogen	Lactofen with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	LTF
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Lactofen is a dark brown to tan solid. Insoluble in water. Used as an herbicide.
Keywords	LACTOFEN; 77501-63-4; 2-Ethoxy-1-methyl-2-oxoethyl 5-(2-chloro-4-trifluoromethylphenoxy)-2-nitrobenzoate; (1-ethoxy-1-oxopropan-2-yl) 5-[2-chloro-4-(trifluoromethyl)phenoxy]-2-nitrobenzoate; Cobra

GENE INFORMATION

Synonyms	LACTOFEN; 77501-63-4; 2-Ethoxy-1-methyl-2-oxoethyl 5-(2-chloro-4-trifluoromethylphenoxy)-2-nitrobenzoate; (1-ethoxy-1-oxopropan-2-yl) 5-[2-chloro-4-(trifluoromethyl)phenoxy]-2-nitrobenzoate; Cobra
----------	--