



Mouse Anti-Fenoxycarb monoclonal antibody, clone FOC (CABT-ZB071)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Fenoxycarb
Immunogen	Fenoxycarb with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	FOC
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Fenoxycarb is a synthetic carbamate ester compound and juvenile hormone mimic that is used as a pesticide. It is characterized as a colorless to white crystalline solid, and exposure occurs by inhalation, ingestion, or contact.
Keywords	FENOXYCARB; 72490-01-8; Insegar; Ethyl (2-(4-phenoxyphenoxy)ethyl)carbamate; ethyl [2-(4-phenoxyphenoxy)ethyl]carbamate

GENE INFORMATION

Synonyms	FENOXYCARB; 72490-01-8; Insegar; Ethyl (2-(4-phenoxyphenoxy)ethyl)carbamate; ethyl [2-(4-phenoxyphenoxy)ethyl]carbamate
----------	---