



# Mouse Anti-Metoprolol monoclonal antibody, clone MPL (CABT-ZB073)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Metoprolol
Immunogen	Metoprolol with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	MPL
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

# BACKGROUND

Introduction	Metoprolol is a cardioselective beta-blocker that is widely used in the treatment of hypertension and angina pectoris. Metoprolol has been linked to rare cases of drug induced liver injury.
Keywords	metoprolol; 51384-51-1; (RS)-Metoprolol; Beatrolol; Spesicor

# GENE INFORMATION

Synonyms	metoprolol; 51384-51-1; (RS)-Metoprolol; Beatrolol; Spesicor
----------	--