



Mouse Anti-Repaglinide monoclonal antibody, clone RGL (CABT-ZB074)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Repaglinide
Immunogen	Repaglinide with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	RGL
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction Repaglinide is a benzoic acid derivative that stimulates insulin secretion from the pancreas and is used in the therapy of type 2 diabetes. Repaglinide has been linked to rare instances of clinically apparent acute liver injury.

Keywords Repaglinide; 135062-02-1; Prandin; NovoNorm; GlucoNorm

GENE INFORMATION

Synonyms Repaglinide; 135062-02-1; Prandin; NovoNorm; GlucoNorm
