



Mouse Anti-Salicylic acid monoclonal antibody, clone SA (CABT-ZB084)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Salicylic acid
Immunogen	Salicylic acid with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	SA
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Salicylic Acid is a beta hydroxy acid that occurs as a natural compound in plants. It has direct activity as an anti-inflammatory agent and acts as a topical antibacterial agent due to its ability to promote exfoliation.

Keywords

Salicylic Acid; $C_6H_4(OH)COOH$; Aspirin

GENE INFORMATION

Synonyms

Salicylic Acid; $C_6H_4(OH)COOH$; Aspirin
