



Mouse Anti-Horse CCL5 monoclonal antibody, clone 172 (CABT-ZB094)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	CCL5
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Horse
Clone	172
Conjugate	Unconjugated
Applications	ELISA Evaluating innate immune responses in horses
Buffer	PBS
Preservative	Sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Chemokine (C-C motif) ligand 5 (also CCL5) is a protein which in humans is encoded by the CCL5 gene. It is also known as RANTES (regulated on activation, normal T cell expressed and secreted).
Keywords	CCL5; chemokine (C-C motif) ligand 5; SISd; eoCP; SCYA5; RANTES; TCP228; D17S136E;

GENE INFORMATION

Synonyms

CCL5; chemokine (C-C motif) ligand 5; SISd; eoCP; SCYA5; RANTES; TCP228; D17S136E;
SIS-delta; C-C motif chemokine 5
