



Mouse Anti-Horse IL-17A monoclonal antibody, clone 76 (CABT-ZB098)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	IL-17A
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Horse
Clone	76
Conjugate	Unconjugated
Applications	ELISA, FACS
Buffer	PBS
Preservative	Sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Interleukin-17A is a protein that in humans is encoded by the IL17A gene. In rodents, IL-17A used to be referred to as CTLA8, after the similarity with a viral gene (O40633).
Keywords	Interleukin 17A; IL17; CTLA8; IL-17; ILA17; CTLA-8; IL-17A

GENE INFORMATION

Synonyms

Interleukin 17A; IL17; CTLA8; IL-17; ILA17; CTLA-8; IL-17A
