



# Mouse Anti-Horse IgA monoclonal antibody, clone 135 (CABT-ZB102)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	IgA
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Horse
Clone	135
Conjugate	Unconjugated
Applications	WB, FACS
Buffer	PBS
Preservative	Sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

## BACKGROUND

Introduction	Immunoglobulin A is an antibody that plays a role in the immune function of mucous membranes. The amount of IgA produced in association with mucosal membranes is greater than all other types of antibody combined.
Keywords	IgA; Immunoglobulin A; Immunoglobulin

# GENE INFORMATION

## Synonyms

IgA; Immunoglobulin A; Immunoglobulin

---