



Rabbit Anti-Human CIB2 monoclonal antibody, clone S112 (CABT-ZB1033)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human CIB2.
Target	CIB2
Immunogen	Recombinant Human CIB2 protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S112
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, ELISA(det), IHC-P We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB706 - CABT-ZB1033 This antibody will detect CIB2 in antibody pair set. [ABPR-ZB286]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CIB2.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Calcium and integrin-binding protein 2 (CIB2) belongs to a protein family with four known members, CIB1 through CIB4, which are characterized by multiple calcium-binding EF-hand domains. Sensorineural hearing loss is genetically heterogeneous. The mutations in CIB2, which encodes a calcium- and integrin-binding protein, are associated with nonsyndromic deafness (DFNB48) and Usher syndrome type 1J (USH1J). Furthermore, in zebrafish and <i>Drosophila melanogaster</i> , CIB2 is essential for the function and proper development of hair cells and retinal photoreceptor cells. We also show that CIB2 is a new member of the vertebrate Usher interactome. Variants in CIB2 can underlie either Usher syndrome type I (USH1J) or nonsyndromic hearing impairment (NSHI) (DFNB48). CIB2 is widely expressed in various human and animal tissues, mainly in skeletal muscle, nervous tissue, inner ear, and retina. The CIB2 protein is responsible for maintaining Ca(2+) homeostasis in cells and interacting with integrins-transmembrane receptors essential for cell adhesion, migration, and activation of signaling pathways. Calcium signaling pathway is crucial for signal transduction in the inner ear, and integrins regulate hair cell differentiation and maturation of the stereocilia.
Keywords	CIB2; calcium and integrin binding family member 2; KIP2; USH1J

GENE INFORMATION

Synonyms	CIB2; calcium and integrin binding family member 2; KIP2; USH1J; DFNB48; calcium and integrin-binding family member 2; KIP 2; 2810434I23Rik; Usher syndrome 1J (autosomal recessive); DNA-dependent protein kinase catalytic subunit-interacting protein 2
Entrez Gene ID	10518
UniProt ID	O75838