



Mouse Anti-Human IL-29 monoclonal antibody, clone NN14 (CABT-ZB1034)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human IL-29.
Target	IL29
Immunogen	Recombinant Human IL-29 Protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	NN14
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB707 - CABT-ZB1034 This antibody will detect IL-29 in antibody pair set. [ABPR-ZB287]
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human IL-29. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Format	Purified, Liquid
Concentration	Lot specific

Size	50 µL, 100 µL, 200 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Interleukin-29 (IL-29), also known as cytokine Zcyto21, Interferon lambda-1, IFN-lambda-1, and IFNL1, is a secreted protein that belongs to the IL-28/IL-29 family. IL-29 is a cytokine with immunomodulatory activity. IL-29 is highly similar in amino acid sequence to the IL-28. IL-28 and IL-29 are induced by a viral infection and showed antiviral activity. IL-28 and IL-29 interacted with a heterodimeric class II cytokine receptor that consisted of IL-10 receptor beta (IL-10Rbeta) and an orphan class II receptor chain, designated IL-28Ralpha. IL-29 plays an important role in host defenses against microbes and its gene is highly upregulated in cells infected with viruses. IL-29 may play a role in antiviral immunity. IL-29 up-regulates MHC class I antigen expression. It is a Ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IL28RA. The ligand/receptor complex seems to signal through the Jak-STAT pathway.
Keywords	IL2RG; interleukin 2 receptor, gamma; P64; CIDX

GENE INFORMATION

Synonyms	IL2RG; interleukin 2 receptor, gamma; P64; CIDX; IMD4; CD132; SCIDX; IL-2RG; SCIDX1; cytokine receptor common subunit gamma; gammaC; CD132 antigen; IL-2R subunit gamma; IL-2 receptor subunit gamma; common cytokine receptor gamma chain
Entrez Gene ID	16186
UniProt ID	P34902