



# Mouse Anti-EIV H3 hemagglutinin monoclonal antibody, clone D-3713 (CABT-ZB105)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	IAV
Immunogen	The hemagglutinin of the prototype equine-2 virus, A/equine/Miami/1/63 (H3N8).
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	EIV
Clone	D-3713
Conjugate	Unconjugated
Applications	IA
Buffer	PBS
Preservative	Sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

## BACKGROUND

Introduction	Influenza A virus (IAV) causes influenza in birds and some mammals, and is the only species of the genus Alphainfluenzavirus of the virus family Orthomyxoviridae. Strains of all subtypes of influenza A virus have been isolated from wild birds, although disease is uncommon. Some isolates of influenza A virus cause severe disease both in domestic poultry and, rarely, in
--------------	--

humans.

---

**Keywords**

HA; Hemagglutinin; Hemagglutinin HA; IAV Influenza A

---

## GENE INFORMATION

**Synonyms**

HA; Hemagglutinin; Hemagglutinin HA; IAV Influenza A

---