



Rabbit Anti-Mouse MCP-1/CCL2 monoclonal antibody, clone S112 (CABT-ZB1050)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Mouse MCP-1/CCL2
Target	CCL2
Immunogen	Recombinant Mouse MCP-1/CCL2 Protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	S112
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB731 - CABT-ZB1050 This antibody will detect MCP-1/CCL2 in antibody pair set. [ABPR-ZB311]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse MCP-1/CCL2.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Monocyte chemoattractant protein 1 (MCP-1), also called CCL2, belongs to a group of CC chemokines located in chromosome 17q11.2. MCP-1 protein interacts with chemokine C-C motif receptor 2 (CCR2) to activate and recruit monocytes, macrophages, CD4+ T cells and immature dendritic cells to the site of infection. The presence of MCP-1 protein in an adequate concentration is important for granuloma formation and M. tuberculosis clearance.
Keywords	CCL2; chemokine (C-C motif) ligand 2; HC11; MCAF

GENE INFORMATION

Synonyms	CCL2; chemokine (C-C motif) ligand 2; HC11; MCAF; MCP1; MCP-1; SCYA2; GDCF-2; SMC-CF; HSMCR30; C-C motif chemokine 2; small-inducible cytokine A2; monocyte secretory protein JE; monocyte chemotactic protein 1; monocyte chemoattractant protein 1; monocyte chemoattractant protein-1; monocyte chemotactic and activating factor; small inducible cytokine subfamily A (Cys-Cys); member 2; small inducible cytokine A2 (monocyte chemotactic protein 1; homologous to mouse Sig-je); anti-MCP-1
Entrez Gene ID	20296
UniProt ID	P10148