



Rabbit Anti-Human ADAM15 monoclonal antibody, clone S118 (CABT-ZB1062)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human ADAM15
Target	ADAM15
Immunogen	Recombinant Human ADAM15 protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S118
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, ELISA(det), IHC-P, FC, ICC/IF We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB745 - CABT-ZB1062 This antibody will detect ADAM15 in antibody pair set. [ABPR-ZB325]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human ADAM15.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	ADAM15, also known as Metargidin, is a type I transmembrane glycoprotein belonging to the ADAM (A Disintegrin and Metalloprotease Domain) family of proteins and is widely expressed in different tissues and cell types. Members of this family contain an amino-terminal metalloprotease domain followed by a disintegrin domain, a cysteine-rich region and a membrane proximal EGF-like domain. The disintegrin domain of ADAM15/metargidin contains an RGD tripeptide sequence, suggesting that it may potentially interact with the integrin family of proteins. ADAM15 is a transmembrane multi-domain proteins implicated in proteolysis, cell-cell and cell-matrix interactions in various disease conditions. There is also evidence supporting a role for ADAM15 in angiogenesis and angioinvasion of tumor cells, which are critical for unrestrained tumor growth and metastatic spread. Given its diverse functions, ADAM15 may represent a pivotal regulatory component of tumor progression, an important target for therapeutic intervention, or emerge as a biomarker of disease progression.
---------------------	---

Keywords	MDC15
-----------------	-------

GENE INFORMATION

Synonyms	MDC15
Entrez Gene ID	8751
UniProt ID	Q13444