



Rabbit Anti-Human CD147 monoclonal antibody, clone S229 (CABT-ZB1076)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human CD147
Target	BSG
Immunogen	Recombinant Human CD147/EMMPRIN/Basigin Protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S229
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB768 - CABT-ZB1076 This antibody will detect CD147 in antibody pair set. [ABPR-ZB348]
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD147 / EMMPRIN / Basigin.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL, 1 mL

Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction CD147/EMMPRIN (Extracellular Matrix Metalloproteinase Inducer), also known as Basigin (BSG), is a transmembrane glycoprotein with different forms resulted from different modes of glycosylation and N-terminal sequence variants. It is a member of the immunoglobulin superfamily with homology to both the immunoglobulin V domain and MHC class II antigen beta-chain. This protein play important roles in variety of events including spermatogenesis, embryo implantation, neural network formation. CD147 induces the production and release of matrix metalloproteinases (MMP) in the surrounding mesenchymal cells and tumor cells, and thereby promotes invasion, metastasis, growth and survival of malignant cells. Furthermore, CD147 also serves as a receptor for extracellular cyclophilinthe and its association with integrins might be important in signal transduction. Recently, CD147 displays increased expression in many cancers, and it has been previously demonstrated to participate in cancer metastasis and progression. Thus, CD147 and its antibody are used as an effective treatment for malignant cancers.

Keywords BSG; basigin; HT-7; CD147

GENE INFORMATION

Synonyms BSG; basigin; HT-7; CD147; EMMPRIN; AI115436; AI325119; gp 42; HT7 antigen; neurothelin; 5A11/Basigin; membrane glycoprotein gp42; basic immunoglobulin superfamily

Entrez Gene ID [682](#)

UniProt ID [P35613](#)