



Mouse Anti-Human Cystatin S monoclonal antibody, clone NN21 (CABT-ZB1085)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human Cystatin S
Target	CST4
Immunogen	Recombinant Human Cystatin S/CST4 Protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	NN21
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA(det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB780 - CABT-ZB1085 This antibody will detect Cystatin S in antibody pair set. [ABPR-ZB360]
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human Cystatin S / CST4. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Format	Purified, Liquid
Concentration	Lot specific

Size	50 µL, 100 µL, 200 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Cystatin-S, also known as Cystatin-4, Salivary acidic protein 1, Cystatin-SA-III and CST4, is a secreted protein which belongs to the cystatin family. Cystatin-4/CST4 is expressed in submandibular and sublingual saliva but not in parotid saliva (at protein level). It is also expressed in saliva, tears, urine and seminal fluid. The cystatin superfamily encompasses proteins that contain multiple cystatin-like sequences. Some of the members are active cysteine protease inhibitors, while others have lost or perhaps never acquired this inhibitory activity. There are three inhibitory families in the superfamily, including the type 1 cystatins (stefins), type 2 cystatins and the kininogens. The type 2 cystatin proteins are a class of cysteine proteinase inhibitors found in a variety of human fluids and secretions. Cystatin-4/CST4 strongly inhibits papain and ficin, partially inhibits stem bromelain and bovine cathepsin C, but does not inhibit porcine cathepsin B or clostripain. Papain is inhibited non-competitively. Cystatin-4/CST4 is an S-type cystatin, based on its high level of expression in saliva, tears and seminal plasma. The specific role in these fluids is unclear but antibacterial and antiviral activity is present, consistent with a protective function.
Keywords	CST4; cystatin S; cystatin-S; cystatin 4

GENE INFORMATION

Synonyms	CST4; cystatin S; cystatin-S; cystatin 4; cystatin-4; cystatin-SA-III; salivary acidic protein 1
Entrez Gene ID	1472
UniProt ID	P01036