



Mouse Anti-Human B7-H3 monoclonal antibody, clone NN17 (CABT-ZB1089)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It reacts with Human B7-H3 It has no cross-reactivity in ELISA with Human cell lysate (293 cell line).
Target	CD276
Immunogen	Recombinant Human B7-H3/CD276 protein
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	NN17
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, ELISA(det), FC, ICC/IF We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB786 - CABT-ZB1089 This antibody will detect B7-H3 in antibody pair set. [ABPR-ZB366]
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human B7-H3 / CD276. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Format	Purified, Liquid

Concentration	Lot specific
Size	50 µL, 100 µL, 200 µL, 1 mL
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction B7-H3 is a member of the B7 family of immune regulatory ligands that is thought to attenuate peripheral immune responses through co-inhibition. It plays an important role in adaptive immune responses, and was shown to either promote or inhibit T-cell responses in various experimental systems. B7-H3 may play an important role in muscle-immune interactions, providing further evidence of the active role of muscle cells in local immunoregulatory processes. B7-H3 is a novel protein structurally related to the B7 family of ligands by the presence of a single set of immunoglobulin-V-like and immunoglobulin-C-like (VC) domains. Previous studies have correlated its overexpression with poor prognosis and decreased tumor-infiltrating lymphocytes in various carcinomas including uterine endometrioid carcinomas, and mounting evidence supports an immuno-inhibitory role in ovarian cancer prognosis. Recently, B7-H3 expression has been reported in several human cancers indicating an additional function of B7-H3 as a regulator of antitumor immunity.

Keywords CD276; CD276 antigen; B7h3; B7RP-2

GENE INFORMATION

Synonyms CD276; CD276 antigen; B7h3; B7RP-2; AU016588; 6030411F23Rik; B7-H3; B7 homolog 3; costimulatory molecule

Entrez Gene ID [80381](#)

UniProt ID [Q5ZPR3](#)