



# Rabbit Anti-Human IL-12 monoclonal antibody, clone S181 (CABT-ZB1092)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |                                                                                                                                                                                      |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Specificity</b>        | It reacts with Human IL-12                                                                                                                                                           |
| <b>Target</b>             | IL12                                                                                                                                                                                 |
| <b>Immunogen</b>          | Recombinant Human IL12/IL12A & IL12B Protein                                                                                                                                         |
| <b>Isotype</b>            | IgG                                                                                                                                                                                  |
| <b>Source/Host</b>        | Rabbit                                                                                                                                                                               |
| <b>Species Reactivity</b> | Human                                                                                                                                                                                |
| <b>Clone</b>              | S181                                                                                                                                                                                 |
| <b>Purification</b>       | Protein A purified                                                                                                                                                                   |
| <b>Conjugate</b>          | Unconjugated                                                                                                                                                                         |
| <b>Applications</b>       | ELISA(det)<br>We recommend the following for sandwich ELISA (Capture - Detection):<br>CABT-ZB790 - CABT-ZB1092<br>This antibody will detect IL-12 in antibody pair set. [ABPR-ZB370] |
| <b>Preparation</b>        | This antibody was obtained from a rabbit immunized with purified, recombinant Human IL12 / IL12A & IL12B.                                                                            |
| <b>Format</b>             | Purified, Liquid                                                                                                                                                                     |
| <b>Concentration</b>      | Lot specific                                                                                                                                                                         |
| <b>Size</b>               | 50 µL, 100 µL, 1 mL                                                                                                                                                                  |

|                     |                                                                                                                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Buffer</b>       | PBS                                                                                                                                                                                                                                                |
| <b>Preservative</b> | None                                                                                                                                                                                                                                               |
| <b>Storage</b>      | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |
| <b>Ship</b>         | Wet ice                                                                                                                                                                                                                                            |

## BACKGROUND

**Keywords** IL-12b; Interleukin 12b; p40; IL-12b

## GENE INFORMATION

**Synonyms** IL-12b; Interleukin 12b; p40; IL-12b; IL12p40; IL-12p40

**Entrez Gene ID** [3593](#)

**UniProt ID** [P29460](#)