



Rabbit Anti-Human RNGTT polyclonal antibody (CABT-ZB110)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross-Reactivity: Human, Mouse, Rat
Target	RNGTT
Immunogen	Recombinant full length Human RNGTT.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC, WB Recommended Dilution: Western Blot: 1:500-1:2000 Immunohistochemistry: 1:50-1:200 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µL
Buffer	PBS with 50% Glycerol

Preservative	0.1% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	RNGTT; RNA guanylyltransferase and 5-phosphatase; HCE; HCE1; hCAP; CAP1A; mRNA-capping enzyme; HCAP1
-----------------	--

GENE INFORMATION

Synonyms	RNGTT; RNA guanylyltransferase and 5-phosphatase; HCE; HCE1; hCAP; CAP1A; mRNA-capping enzyme; HCAP1
-----------------	--

Entrez Gene ID	8732
-----------------------	----------------------

UniProt ID	O60942
-------------------	------------------------