



# Rabbit Anti-Human RNGTT polyclonal antibody (CABT-ZB111)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	RNGTT
Immunogen	Synthetic peptide corresponding to the N-terminal region of Rngtt.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat, Bovine, Dog, Guinea pig, Horse, Rabbit
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µL
Buffer	PBS with 2% sucrose
Preservative	0.09% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

## BACKGROUND

**Keywords**

RNGTT; RNA guanylyltransferase and 5-phosphatase; HCE; HCE1; hCAP; CAP1A; mRNA-capping enzyme; HCAP1

---

## GENE INFORMATION

**Synonyms**

RNGTT; RNA guanylyltransferase and 5-phosphatase; HCE; HCE1; hCAP; CAP1A; mRNA-capping enzyme; HCAP1

---

**Entrez Gene ID**

[8732](#)

---

**UniProt ID**

[O60942](#)

---