



Rabbit Anti-Human ZDHHC1 (aa 14-43) polyclonal antibody (CABT-ZB114)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	ZDHHC1
Immunogen	A KLH conjugated synthetic peptide between 14-43 amino acids from the N-terminal region of human ZDHHC1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IF, FC Recommended Dilution: WB: 1:1000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µL
Buffer	PBS
Preservative	0.09% sodium azide

Storage	Long time storage is recommended at -20°C.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Keywords	ZDHHC1; zinc finger, DHHC-type containing 1; C16orf1, chromosome 16 open reading frame 1; probable palmitoyltransferase ZDHHC1; HSU90653; ZNF377; DHHC-1; zinc finger protein 377; zinc finger, DHHC domain containing 1; DHHC-domain-containing cysteine-rich protein; zinc finger DHHC domain-containing protein 1; DHHC domain-containing cysteine-rich protein 1; C16orf1
-----------------	---

GENE INFORMATION

Synonyms	ZDHHC1; zinc finger, DHHC-type containing 1; C16orf1, chromosome 16 open reading frame 1; probable palmitoyltransferase ZDHHC1; HSU90653; ZNF377; DHHC-1; zinc finger protein 377; zinc finger, DHHC domain containing 1; DHHC-domain-containing cysteine-rich protein; zinc finger DHHC domain-containing protein 1; DHHC domain-containing cysteine-rich protein 1; C16orf1
-----------------	---

Entrez Gene ID	29800
-----------------------	-----------------------

UniProt ID	Q8WTX9
-------------------	------------------------