



Rabbit Anti-Human ZDHC12 polyclonal antibody (CABT-ZB116)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of ZDHC12 at Human, Mouse, Rat
Target	ZDHC12
Immunogen	Synthesized peptide derived from human ZDHC12
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB Recommended Dilution: WB: 1:500-2000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS with 50% Glycerol

Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	ZDHHC12; zinc finger, DHHC-type containing 12; ZNF400; probable palmitoyltransferase ZDHHC12; DHHC-12; zinc finger protein 400; zinc finger, DHHC domain containing 12; zinc finger DHHC domain-containing protein 12
-----------------	---

GENE INFORMATION

Synonyms	ZDHHC12; zinc finger, DHHC-type containing 12; ZNF400; probable palmitoyltransferase ZDHHC12; DHHC-12; zinc finger protein 400; zinc finger, DHHC domain containing 12; zinc finger DHHC domain-containing protein 12
-----------------	---

Entrez Gene ID	84885
-----------------------	-----------------------

UniProt ID	Q96GR4
-------------------	------------------------