



Rabbit Anti-Human ZDHC11 polyclonal antibody (CABT-ZB117)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of ZDHC11 at Human, Mouse
Target	ZDHC11
Immunogen	Synthesized peptide derived from human ZDHC11
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB Recommended Dilution: WB: 1:500-2000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS with 50% Glycerol

Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	ZDHHC11; zinc finger, DHHC-type containing 11; probable palmitoyltransferase ZDHHC11; FLJ13153; ZNF399; 4933421L13Rik; DHHC 11; DHHC-11; MGC151180; MGC151184
-----------------	---

GENE INFORMATION

Synonyms	ZDHHC11; zinc finger, DHHC-type containing 11; probable palmitoyltransferase ZDHHC11; FLJ13153; ZNF399; 4933421L13Rik; DHHC 11; DHHC-11; MGC151180; MGC151184
-----------------	---

Entrez Gene ID	79844
-----------------------	-----------------------

UniProt ID	Q9H8X9
-------------------	------------------------