



Rabbit Anti-Human ZDHC22 polyclonal antibody (CABT-ZB118)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of ZDHC22 at Human, Mouse, Rat
Target	ZDHC22
Immunogen	Synthesized peptide derived from human ZDHC22
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB Recommended Dilution: WB: 1:500-2000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS with 50% Glycerol

Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	DHHC22; Zinc Finger, DHHC Domain Containing; Putative Palmitoyltransferase ZDHHC22; ZDHHC22
-----------------	---

GENE INFORMATION

Synonyms	DHHC22; Zinc Finger, DHHC Domain Containing; Putative Palmitoyltransferase ZDHHC22; ZDHHC22
-----------------	---

Entrez Gene ID	283576
-----------------------	------------------------

UniProt ID	Q8N966
-------------------	------------------------