



Mouse Anti-GD2 antibody 14G2a monoclonal antibody, clone 1A7 (CABT-ZB1182)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A mouse anti-GD2 antibody 14G2a.
Isotype	IgG1, κ
Source/Host	Mouse
Clone	1A7
Purification	Protein A purified
Conjugate	Unconjugated
Applications	We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB1177 - CABT-ZB1180
Format	Liquid
Concentration	Lot specific
Size	200 µg
Buffer	PBS
Preservative	0.02% Proclin 300
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Ship	Wet ice