



Mouse Anti-D3H44 monoclonal antibody, clone 6A6 (CABT-ZB1184)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Humanised monoclonal antibody D3H44.
Isotype	IgG1, λ
Source/Host	Mouse
Species Reactivity	Human
Clone	6A6
Purification	Protein A purified
Conjugate	Unconjugated
Format	Liquid
Concentration	Lot specific
Size	200 μ g
Buffer	PBS
Preservative	0.02% Proclin 300
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Ship	Wet ice